



IOS CONNECTION GUIDE

Sirona (CEREC)

Go to Customer.ConnectCaseCenter.com

Create Account or Log In to your Connect Case Center account.

Select "My Connect Case Center." Then "Partner."

Enter "Image Dental Arts" in the box for Company Name.

Click FIND.

Check the box to the bottom left for "Image Dental Arts".

Click ADD to complete the registration.

DS Core (Sirona & CEREC will become DS Core on Nov. 14th)

Call 888-268-9095 to register.

In the laboratories section, under preferred labs, search for Image Dental Arts, and request Image Dental Arts be added as a preferred lab.

DS Core will then confirm with Image Dental Arts and add connection remotely.

iTero®

Call 800-577-8767.

Select Option 2.

Request that Image Dental Arts (ID: 3710) is added to your scanner, and identify our lab using our phone number: 828-225-3812

After Image Dental Arts has been added, restart your scanner.

After connecting with us as a lab, select Image Dental Arts on your scanner when sending files.

TRIOS®

Go to us.3shapecommunicate.com in a web browser

Connect with Image Dental Arts as lab by searching imagedigital.ida@gmail.com

After connecting to us as a Lab, select Image Dental Arts when sending files.

Medit Link

Go to the Medit link website and create an account. Log in to your account.

Click on "Partners" on the left side. Type imagedigital.ida@gmail.com in the search bar.

Find Image Dental Arts and click "Request Partnership."

Check the details and press "Request Partnership" again.

We will approve your request quickly, and your status will update to "Partners."

DEXIS IS Connect

Go to <https://dentalconnect.dexis.com> on a web browser and login

Login to DEXIS IS Connect as a Dentist.

In the toolbar, click .

In the Partners window, click Invite a partner.

In the Email field, enter the email address of the Laboratory.

In the Message box, either leave the default text or write your own message.

Click Send.

The invitation message is sent to the Laboratory and a Partner entry is created in the list of Partners.

CARESTREAM CS Connect

Go to <https://csdentalconnect.com/> on a web browser and login

Login to CS Connect as a Dentist.

In the toolbar, click .

In the Partners window, click Invite a partner.

In the Email field, enter the email address of the Laboratory.

In the Message box, either leave the default text or write your own message.

Click Send.

The invitation message is sent to the Laboratory and a Partner entry is created in the list of Partners.

DDX

Go to <https://www.ddxdental.com/> and login

Click the-search-button, this will attempt to help you find the Lab on the network.

Enter the Country (US) and Lab Name (Image Dental Arts). Alternately you can search by State (NC) or

Zip Code (28803) and click the associated the-search-button.

If the Partner Lab you are looking for is in the Search Results, you will need to Request and Account with that Lab. Click the-go-to-lab-button. The Request an Account Page will open, and now you are now on

your Partner Labs site. Click Confirm Preferences, then Click the Request Account button. The Partner Lab's site will now reflect your current relationship with them, which is Account Activation

Shining 3D DentalCloud

Go to <https://dental3dcloud.com> in a web browser

Go to the Connection.

Click New in the top right corner to open the New relationship network dialog.

Enter the institution name (Image Dental Arts) to search for the designated lab. Click +, and a confirmation dialog will appear.

Click Confirm to invite the lab to establish the connection.